**ENGLISH****S-2444™**

For Laboratory Use Only

S-2444™

S-2444 is a chromogenic substrate for urokinase.

COMPOSITION

Each vial contains chromogenic substrate S-2444 25 mg and mannitol 40 mg added as a bulking agent.

CHEMISTRY

Chemical name: L-Pyroglutamyl-glycyl-L-arginine-p-Nitroaniline hydrochloride.

Formula: <Glu-Gly-Arg-pNA · HCl

Mol. wt.: 498.9

$\epsilon_{316\text{ nm}^1}$ 1.27 · 10⁴ mol⁻¹ · L · cm⁻¹

Solubility: > 10 mmol/L in H₂O

Stability: Substance: Stable until expiry date if stored at 2-8°C. Avoid exposure to light. The substance is hygroscopic and should be stored dry. Solution: 2 mmol/L in H₂O is stable for at least 6 months at 2-8°C. Contamination by micro-organisms may cause hydrolysis.

Suitable stock solution:

2-3 mmol/L in H₂O.**PRINCIPLE**

<Glu-Gly-Arg-pNA $\xrightarrow{\text{Enzyme}}$ <Glu-Gly-Arg-OH+pNA

The method for the determination of activity is based on the difference in absorbance (optical density) between the pNA formed and the original substrate. The rate of pNA formation, i.e. the increase in absorbance per second at 405 nm, is proportional to the enzymatic activity and is conveniently determined with a photometer.

KINETIC DATA

Urokinase: Ploug U

$K_m = 9 \cdot 10^{-5}$ mol/L,
 $V = 3.1 \cdot 10^{-10}$ mol/min

CTA U

$K_m = 6 \cdot 10^{-5}$ mol/L,
 $V = 1.3 \cdot 10^{-10}$ mol/min

Determined at 37°C in 2.5 mL of 0.05 mol/L Tris buffer pH 8.8, I 0.05.

STANDARDIZATION

An activity of $\Delta A/\text{min} = 0.05$ (37°C) is obtained by using 0.3 mmol/L substrate and:

1. 25 Ploug U/mL of urokinase (Leo, Copenhagen, Denmark)
2. 40 CTA U/mL of urokinase (Abbott, Chicago, USA)

APPLICATIONS

The substrate has been used for the determination of:

1. Urokinase in purified preparations (1,2)
2. Urokinase inhibitors in plasma (3)
3. Urokinase activity in urine (4)

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301983R0

REFERENCES

1. CLAESON G et al.: Methods for the determination of prekallikrein in plasma glandular kallikrein and urokinase. *Haemostasis* 7, 76-78 (1978).
2. Chromogenix AB: Standardization of urokinase with S-2444. Laboratory Instruction.
3. GALLIMORE M & FAREID E: Studies on human inhibitors of plasmin, plasma kallikrein, trypsin, thrombin and urokinase using chromogenic substrate assays. In *Chromogenic Peptide Substrates*. Ed. M F Scully and V V Kakkar. Churchill Livingstone 248-261 (1979).
4. PAAR D and MARUHN D: Spectrometric determination of urokinase in urine after gel filtration, using the chromogenic substrate S-2444. *J Clin Chem Clin Biochem* 18, 557-562 (1980).

CHROMOGENIX

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Instrumentation Laboratory SpA